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
## 1.3.3 Cleavage


**R<sub>3</sub> (TFA)** is a highly-volatile, strong organic acid. R<sub>3</sub> cleaves the PITC-coupled amino acid residue from the amino terminus of the protein, thereby producing the anilinothiazolinone (ATZ) derivative of the amino acid. The remaining protein is left with a new amino terminus which will be coupled with PITC in the next cycle.

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### 1.3.5 PTH Conversion

**R4A (25% TFA in water – aqueous TFA)** The transferred extraction solution (S3 and/or S2B + ATZ amino acid), along with most of the initial delivery of S4B, is evaporated from the conversion flask. An aliquot of R4A is then delivered to the flask to convert the ATZ-AAs to the more stable PTH-AAs. The presence of dithiothreitol in the R4A minimizes oxidation of the ATZ amino acids. In standard sequencer cycles, conversion occurs for approximately nine minutes at 64 °C. These conditions are sufficient to give virtually complete conversion of all amino acids except glycine (about 15% of derivatized glycine remains in the intermediate PTC-glycine form).

**R5 (acetonitrile)** is used to introduce an aliquot of the PTH standard mix to the flask for subsequent HPLC analysis.

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